## **EUCHEMBI**

## THE EUROPEAN CHEMISTRY AND BIOTECHNOLOGY JOURNAL

Eur. Chem. Biotechnol. J. 1: 39-46 (2024)

www.euchembioj.com

# Research article

# Salt-assisted liquid–liquid microextraction for determination of haloperidol in human plasma by LC-MS/MS

Selen Al<sup>1</sup> ORCID 0009-0009-2014-1998 ayseselen.al@ogr.iu.edu.tr

Aykut Kul<sup>1\*</sup> ORCID 0000-0002-1195-0939 aykut.kul@istanbul.edu.tr

Olcay Sagirli<sup>1</sup> ORCID 0000-0002-0434-2150 osagirli@istanbul.edu.tr

# <sup>1</sup>Department of Analytical Chemistry, Faculty of Pharmacy, Istanbul University, 34452 Istanbul, Türkiye

\*Corresponding author: aykut.kul@istanbul.edu.tr; Ph.: +90-212-440-0000; Fax: +90-212-440-0252.

Received: 20.10.2023 Accepted: 06.12.2023 Published: 02.01.2024

Cite this article: Al, S., Kul, A., & Sagirli, O. (2024). Salt-assisted liquid–liquid microextraction for determination of haloperidol in human plasma by LC-MS/MS. *The European chemistry and biotechnology journal*, *1*, 39-46. https://doi.org/10.62063/ecb-16

## Abstract

Haloperidol is an antipsychotic used in the treatment of schizophrenia. Compared to other antipsychotics, it is widely used in developing countries due to its affordable price. Haloperidol has a narrow therapeutic range and variable pharmacokinetics; therefore, therapeutic drug monitoring (TDM) is recommended. For this reason, in this study, an easily applicable, fast, selective, accurate, reliable, and economical LC-MS/MS method was developed for the determination of haloperidol in human plasma for use in TDM and also method was validated according to European Medicines Agency (EMA) Bioanalytical method validation guidelines. In the developed method, analyte and internal standard were extracted from plasma by salt-assisted liquid-liquid microextraction (SALLME) technique and after that injected to the LC system. The limit of quantification of haloperidol was determined as 1 ng/ml. The calibration curve was validated between 1-15 ng/ml, with correlation coefficients >0.99. In addition, the developed method was used to determine drug concentration levels in the plasma of real patients.

Keywords: haloperidol, LC-MS/MS, SALLME, TDM

#### Introduction

Haloperidol is a butyrophenone derivative antipsychotic agent used in the treatment of schizophrenia and related psychoses (Lee et al., 2007). It was synthesized in 1958 and first marketed in the United States in 1967 (Gerace et al., 2012). While there is a wide choice of antipsychotics, some very expensive, in developed countries, older and less expensive drugs such as haloperidol is still widely used in developing countries (Leucht et al., 2008). Although haloperidol is frequently prescribed and often used as a comparator in clinical trials, both acute and long-term use can cause extrapyramidal

side effects such as dystonia, dyskinesia and parkinsonism (Grubor et al., 2020). This drug exerts its therapeutic effect through D2 receptor antagonism, as with other antipsychotic drugs (Tyler et al., 2017). It has a rapid onset of action and an extensive volume of tissue distribution. It is also lipophilic and highly protein bound (Wang et al. 2012). Only 1% of the administered dose of haloperidol is excreted unchanged in the urine because the drug is extensively metabolized in the liver (Gassó et al., 2013). Its metabolism includes conversion to an "inactive metabolite" (50/60%) by glucuronidation, to "reduced haloperidole" (an active metabolite - 23%) by reduction and back oxidation, and to a "pyridium metabolite" (a toxic metabolite 20/30%) by N-dealkylation (Vella-Brincat and Macleod 2004). The CYP2D6 genotype has an important role in the excretion of "haloperidol" and "reduced haloperidol" and may pose a risk of serious adverse effects. In addition, body weight, and smoking may cause changes in pharmacokinetics (Brockmöller et al., 2002). Haloperidol can be used with a variety of medications such as antipsychotics, antidepressants, and benzodiazepines (Shin et al. 2001). Combination antipsychotic therapy is a common practice, but it also has effects such as drug-drug interactions, variability in serum levels, and difficulties in patient adherence to treatment (Lally and MacCabe 2015). Because of this variable pharmacokinetics and narrow therapeutic range, haloperidol therapy requires personalized optimization. For this reason, therapeutic drug monitoring is recommended for haloperidol (Kudo and Ishizaki 1999). The recommended plasma concentration for haloperidol is 1-10 ng/mL, while the alert level is 15 ng/mL (Hiemke et al., 2018). Therapeutic drug monitoring (TDM) helps the clinician reach targeted concentrations and direct dosage to reduce the risk of side effects (Eliasson et al. 2013). TDM is performed by chromatographic and immunoassay methods, but in immunoassays the presence of various other components in the sample can lead to a lack of analytical specificity (Ates et al., 2020). For this reason, reliable, sensitive, high-quality analytical methods by instruments combined with Ultraviolet (UV), Fluorescent (FL), and mass spectrometry (MS) detectors are the main techniques used for TDM (Tuzimski and Petruczynik 2020).

Salt-assisted liquid liquid microextraction (SALLME) is a simple, fast and easily applicable sample preparation technique that uses various salts and water-miscible organic solvents (Gupta et al. 2009). Currently, there is a lack of any SALLME technique in the literature for the determination of haloperidol in human plasma. Therefore, the aim of this study was to develop a rapid, economical and easily applicable LC-MS/MS method for the determination of haloperidol in human plasma using SALLME technique.

#### Materials and methods

#### Materials and reagents

Haloperidol, haloperidol-d4, sodium hydroxide, sodium chloride, and organic solvents were procured from Merck (Darmstadt, Germany). Water was purified by a Purelab option-Q water system from Elga (Elga, Lane End, UK).

Stock solutions of haloperidol were prepared separately in methanol at a 1 mg/mL concentration. A working standard solution was prepared by diluting the stock standard solutions with an appropriate volume of methanol. Calibration standards and quality control samples were prepared by spiking the (internal standard) IS and working standard solution into blank plasma. The concentrations of the calibration standard solutions were 1, 2, 4, 6, 8, 10, 13, and 15 ng/mL, while the concentrations of the quality control samples were 1, 3, 7.5,12 ng/mL. All solutions were kept at -20°C when not in use.

#### Sample preparation

The salt-assisted liquid-liquid microextraction technique was followed for cleaning the plasma samples. The sample preparation procedure is shown in Figure 1.




Figure 1. Sample preparation.

#### LC-MS/MS analysis

The method was performed on an Agilent 1200 series HPLC system coupled to an Agilent 6460 triple quadrupole tandem mass spectrometry with an InterSustain C18 (3.0  $\mu$ m, 3.0× 100 mm) column. The column temperature and the injection volume were 30 °C and 5  $\mu$ L respectively. The isocratic elution was performed using mobile phase A (water with 0.1% formic acid) and mobile phase B (acetonitrile with 0.1% formic acid) (60:40). The total run time was 4 min with flow rate of 0.4 mL/min. The desolvation gas temperature, nebulizer, and desolvation gas flow were set at 300 °C, 45 psi, 11 L/min. The typical ion transitions were 376.2 $\rightarrow$  165.1 for haloperidol and m/z 380.2 $\rightarrow$ 169.0 for haloperidol-d4. Analytes were detected using electrospray ionization in the positive mode.

#### Validation of the method

Method validation was performed to evaluate the selectivity, lower limit of quantification (LLOQ), precision, accuracy, dilution integrity, matrix effects, carryover, calibration curve performance, and stability, according to European Medicines Agency (EMA) Bioanalytical Method Validation Guideline (EMA 2011). Selectivity is the ability of the method to separate the analyte from other substances in the matrix. For this purpose, blank plasmas from six different sources were analyzed and examined for any interference. The lower limit of quantification (LLOQ) is the lowest concentration of analyte that can be measured in a sample with acceptable accuracy, reliability, and precision. LLOQ was determined by considering the therapeutic range of the analyte. Accuracy and precision parameters were assessed by preparing quality control samples at 4 different concentrations (LLOQ, low-QC, medium-QC, and high-QC) and analyzing them six times on three different days. Calibration curve performance was evaluated using 8 different calibration concentration levels for the analyte and the calibration curve was defined between the lower limit of quantification and upper limit of quantification (ULOQ). Carryover was evaluated by injecting a blank sample after the standard at the upper limit of quantification. Dilution integrity was assessed by diluting a sample with a higher concentration than the ULOQ sample to within the range of the calibration curve. The stability parameter was evaluated using freshly prepared samples of the analytes and internal standard under 3 different conditions (4 weeks at -20°C, 24 hours in autosampler. 3 times freeze-thaw). In the matrix effect study, low and high QC samples were prepared in the plasma matrix and the solution obtained from six different sources was analyzed. The matrix effect was calculated by proportioning the response values obtained.

#### Sample analysis

The developed method has been successfully applied to plasma samples from patients on drug therapy.

## **Results and discussion**

#### Method optimization

In this study, various amounts of sodium chloride salt (150, 200, and 300 mg), organic solution types (acetonitrile and ethyl acetate), and volumes (300, 400, and 500  $\mu$ L) were tested for phase separation. 200 mg of sodium chloride and 300  $\mu$ L acetonitrile produced the best yield.

#### Method validation

No interference was detected in the retention time of the analytes and internal standard in blank plasma. The retention times for the internal standard and haloperidol (Figure 2). The lower limit of quantification of haloperidol was determined as 1 ng/mL. The calibration curves were validated between 1-15 ng/mL for haloperidol, with correlation coefficients >0.99. No interference was detected in the blank plasma injected after the ULOQ sample.




Figure 2. Chromatograms of the haloperidol (a) and IS (b).

All accuracy and precision values are given in Table 1. It is seen that the results comply with the EMA criteria.

Table 1. Accuracy and precision data for haloperidol.

QC concentration	Intra-assay (n=6)		Inter-assay (n=18)		
(ng/mL)	Accuracy (%)	CV (%)	Accuracy (%)	CV (%)	
1 (LLOQ)	102.9	4.4	96.3	4.2	
3 (low)	103.1	2.1	106.8	2.9	
7.5 (medium)	106.4	3.8	100.9	3.5	
12 (high)	102.5	5.7	106.3	4.8	

QC: quality control, LLOQ: lower limit of quantification, CV: coefficient of variation

The accuracy and precision of the diluted samples were 101.7% and 3.4% for haloperidol. The accuracy values of the stability analyses for haloperidol and the internal standard are shown in Table 2. The results comply with EMA acceptance criteria.

#### Table 2. The stability of haloperidol in plasma.

		Stability (% Accuracy	/)
QC concentration (ng/mL)	24 hours in the autosampler (n=6)	4 weeks at −20 °C (n=6)	three freeze-thaw cycles (n=6)
3 (low)	104.3	102.9	108.7
12 (high)	105.1	103.1	106.5

QC: quality control

The matrix effect results were within EMA's acceptable ranges (CV, 15%) (Table 3).

Table 3. The matrix effect data of haloperidol.

QC concentration (ng/mL)	*Normalized matrix effect factor (n=6)				CV (%)		
3 (low)	1.0	1.0	1.0	1.1	1.0	1.1	5.0
12 (high)	1.0	1.1	1.0	1.0	0.9	1.0	6.3
Or availter a sustant OV/s as	affiniant of						

QC: quality control, CV: coefficient of variation

#### Application to real sample

Five actual samples were successfully used to test the developed method. Table 4 presents all results.

 Table 4. The results of the samples.

Sample	Age Dose of haloperidol		Concentration of plasma (ng/mL)		
Sample-1	19	5 mg	5.5		
Sample-2	44	10 mg	2.4		
Sample-3	14	10 mg	9.9		
Sample-4	46	10 mg	2.3		
Sample-5	24	10 mg	1.4		

#### Method comparison

There are various LC-MS/MS methods in the literature for the determination of haloperidol in plasma and serum. These methods have various sample preparation techniques such as protein precipitation (Juenke et al. 2013, Gradinaru et al., 2014, Domingues et al., 2016, Cao et al. 2020), liquid-liquid extraction (Remane et al., 2011, Montenarh et al., 2016, Degreef et al., 2021), and solid phase extraction (Hempenius et al., 1999, Khelfi et al., 2018). Among these sample preparation techniques, protein precipitation is relatively easy but it has a low recovery. Liquid-liquid extraction requires the use of a lot of organic solvents, while solid-phase extraction is quite expensive compared to these methods. The Sample preparation technique presented in this study is not only easy to apply but also more environmentally friendly as it uses less organic solvents.

## Conclusions

The LC-MS/MS method developed using the SALLME technique is simple, accurate, and precise. The values of the European Medicines Agency's bioanalytical method validation parameters are acceptable and therefore, the method is suitable for the therapeutic drug monitoring of haloperidol. Additionally, the SALLME technique can be used as an extraction method in many biotechnological studies.

## Acknowledgements

Not applicable.

## Funding

Not applicable.

#### **Conflict of interest**

The authors declare no conflict of interest.

## Data availability statement

Data can be obtained from the corresponding author upon a reasonable request.

#### Ethics committee approval

The study was approved by the Ethics Committee of the University of Üsküdar (No. 7 of 29.07.2021) and was conducted in accordance with the Helsinki Declaration.

#### Authors' contribution statement

The authors acknowledge their contributions to this paper as follows: **Study conception and design**: O.S., A.K.; **Data collection**: S.A.; **Analysis and interpretation of results**: S.A.; **Manuscript draft preparation**: A.K. All authors reviewed the results and approved the final version of the manuscript.

#### References

Agency, E.M., 2011. Guideline on bioanalytical method validation. EMEA/CHMP/EWP/192217/2009.

Ates, H. C., Roberts, J. A., Lipman, J., Cass, A. E. G., Urban, G. A., & Dincer, C. (2020). On-site therapeutic drug monitoring. *Trends in biotechnology*, *38*(11), 1262-1277. https://doi.org/10.1016/j.tibtech.2020.03.001

Brockmöller, J., Kirchheiner, J., Schmider, J., Walter, S., Sachse, C., Müller-Oerlinghausen, B., & Roots, I. (2002). The impact of the CYP2D6 polymorphism on haloperidol pharmacokinetics and on the outcome of haloperidol treatment. *Clinical pharmacology and therapeutics*, 72(4), 438-452. https://doi.org/10.1067/mcp.2002.127494

Cao, Y., Zhao, F., Chen, J., Huang, T., Zeng, J., Wang, L., Sun, X., Miao, Y., Wang, S., & Chen, C. (2020). A simple and rapid LC-MS/MS method for the simultaneous determination of eight antipsychotics in human serum, and its application to therapeutic drug monitoring. Journal of chromatography. B, Analytical technologies in the biomedical and life sciences, 1147, 122129. https://doi.org/10.1016/j.jchromb.2020.122129

Degreef, M., Berry, E. M., Maudens, K. E. K., & van Nuijs, A. L. N. (2021). Multi-analyte LC-MS/MS quantification of 38 antipsychotics and metabolites in plasma: Method validation & application to routine analyses. *Journal of chromatography. B, Analytical technologies in the biomedical and life sciences*, *1179*, 122867. https://doi.org/10.1016/j.jchromb.2021.122867

Domingues, D. S., Pinto, M. A., de Souza, I. D., Hallak, J. E., Crippa, J. A., & Queiroz, M. E. (2016). Determination of Drugs in Plasma Samples by High-Performance Liquid Chromatography-Tandem

Mass Spectrometry for Therapeutic Drug Monitoring of Schizophrenic Patients. *Journal of analytical toxicology*, *40*(1), 28-36. https://doi.org/10.1093/jat/bkv107

Eliasson, E., Lindh, J. D., Malmström, R. E., Beck, O., & Dahl, M. L. (2013). Therapeutic drug monitoring for tomorrow. *European journal of clinical pharmacology*, 69 Suppl 1, 25-32. https://doi.org/10.1007/s00228-013-1504-x

Gassó, P., Papagianni, K., Mas, S., de Bobadilla, R. F., Arnaiz, J. A., Bernardo, M., & Lafuente, A. (2013). Relationship between CYP2D6 genotype and haloperidol pharmacokinetics and extrapyramidal symptoms in healthy volunteers. *Pharmacogenomics*, *14*(13), 1551-1563. https://doi.org/10.2217/pgs.13.150

Gerace, E., Salomone, A., Pellegrino, S., & Vincenti, M. (2012). Evidence of Haldol (haloperidol) longterm intoxication. *Forensic science international*, *215*(1-3), 121-123. https://doi.org/10.1016/j.forsciint.2011.01.038

Gradinaru, J., Vullioud, A., Eap, C. B., & Ansermot, N. (2014). Quantification of typical antipsychotics in human plasma by ultra-high performance liquid chromatography tandem mass spectrometry for therapeutic drug monitoring. *Journal of pharmaceutical and biomedical analysis*, *88*, 36-44. https://doi.org/10.1016/j.jpba.2013.07.041

Grubor, M., Zivkovic, M., Sagud, M., Nikolac Perkovic, M., Mihaljevic-Peles, A., Pivac, N., Muck-Seler, D., & Svob Strac, D. (2020). HTR1A, HTR1B, HTR2A, HTR2C and HTR6 Gene Polymorphisms and Extrapyramidal Side Effects in Haloperidol-Treated Patients with Schizophrenia. *International journal of molecular sciences*, *21*(7), 2345. https://doi.org/10.3390/ijms21072345

Gupta, M., Jain, A., & Verma, K. K. (2009). Salt-assisted liquid-liquid microextraction with water-miscible organic solvents for the determination of carbonyl compounds by high-performance liquid chromatography. *Talanta*, *80*(2), 526-531. https://doi.org/10.1016/j.talanta.2009.07.021

Hempenius, J., Steenvoorden, R. J., Lagerwerf, F. M., Wieling, J., & Jonkman, J. H. (1999). 'High throughput' solid-phase extraction technology and turbo ionspray LC-MS-MS applied to the determination of haloperidol in human plasma. *Journal of pharmaceutical and biomedical analysis*, 20(6), 889-898. https://doi.org/10.1016/s0731-7085(99)00107-7

Hiemke, C., Bergemann, N., Clement, H. W., Conca, A., Deckert, J., Domschke, K., Eckermann, G., Egberts, K., Gerlach, M., Greiner, C., Gründer, G., Haen, E., Havemann-Reinecke, U., Hefner, G., Helmer, R., Janssen, G., Jaquenoud, E., Laux, G., Messer, T., Mössner, R., ... Baumann, P. (2018). Consensus Guidelines for Therapeutic Drug Monitoring in Neuropsychopharmacology: Update 2017. *Pharmacopsychiatry*, *51*(1-02), e1. https://doi.org/10.1055/s-0037-1600991

Juenke, J. M., Brown, P. I., Urry, F. M., Johnson-Davis, K. L., & McMillin, G. A. (2013). Simultaneous UPLC-MS/MS assay for the detection of the traditional antipsychotics haloperidol, fluphenazine, perphenazine, and thiothixene in serum and plasma. *Clinica chimica acta; international journal of clinical chemistry*, *423*, 32-34. https://doi.org/10.1016/j.cca.2013.04.014

Khelfi, A., Azzouz, M., Abtroun, R., Reggabi, M., & Alamir, B. (2018). Determination of Chlorpromazine, Haloperidol, Levomepromazine, Olanzapine, Risperidone, and Sulpiride in Human Plasma by Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS). *International journal of analytical chemistry*, 2018, 5807218. https://doi.org/10.1155/2018/5807218

Kudo, S., & Ishizaki, T. (1999). Pharmacokinetics of haloperidol: an update. *Clinical pharmacokinetics*, *37*(6), 435-456. https://doi.org/10.2165/00003088-199937060-00001

Lally, J., & MacCabe, J. H. (2015). Antipsychotic medication in schizophrenia: a review. *British medical bulletin*, *114*(1), 169-179. https://doi.org/10.1093/bmb/ldv017

Lee, Y., Wang, P. K., Lai, H. Y., Yang, Y. L., Chu, C. C., & Wang, J. J. (2007). Haloperidol is as effective as ondansetron for preventing postoperative nausea and vomiting. *Canadian journal of anaesthesia = Journal canadien d'anesthesie*, *54*(5), 349–354. https://doi.org/10.1007/BF03022656

Leucht, C., Kitzmantel, M., Chua, L., Kane, J., & Leucht, S. (2008). Haloperidol versus chlorpromazine for schizophrenia. *The Cochrane database of systematic reviews*, (1), CD004278. https://doi.org/10.1002/14651858.CD004278.pub2

Montenarh, D., Hopf, M., Maurer, H. H., Schmidt, P., & Ewald, A. H. (2016). Development and validation of a multi-analyte LC-MS/MS approach for quantification of neuroleptics in whole blood, plasma, and serum. *Drug testing and analysis*, *8*(10), 1080-1089. https://doi.org/10.1002/dta.1923

Remane, D., Meyer, M. R., Wissenbach, D. K., & Maurer, H. H. (2011). Ultra high performance liquid chromatographic-tandem mass spectrometric multi-analyte procedure for target screening and quantification in human blood plasma: validation and application for 31 neuroleptics, 28 benzodiazepines, and Z-drugs. *Analytical and bioanalytical chemistry*, *401*(4), 1341-1352. https://doi.org/10.1007/s00216-011-5187-9

Shin, J. G., Kane, K., & Flockhart, D. A. (2001). Potent inhibition of CYP2D6 by haloperidol metabolites: stereoselective inhibition by reduced haloperidol. *British journal of clinical pharmacology*, *51*(1), 45-52. https://doi.org/10.1046/j.1365-2125.2001.01313.x

Tuzimski, T., & Petruczynik, A. (2020). Review of Chromatographic Methods Coupled with Modern Detection Techniques Applied in the Therapeutic Drugs Monitoring (TDM). *Molecules (Basel, Switzerland)*, *25*(17), 4026. https://doi.org/10.3390/molecules25174026

Tyler, M. W., Zaldivar-Diez, J., & Haggarty, S. J. (2017). Classics in Chemical Neuroscience:Haloperidol.ACSchemicalneuroscience,8(3),444–453.https://doi.org/10.1021/acschemneuro.7b00018

Vella-Brincat, J., & Macleod, A. D. (2004). Haloperidol in palliative care. *Palliative medicine*, 18(3), 195–201. https://doi.org/10.1191/0269216304pm881oa

Wang, E. H., Mabasa, V. H., Loh, G. W., & Ensom, M. H. (2012). Haloperidol dosing strategies in the treatment of delirium in the critically ill. *Neurocritical care*, *16*(1), 170–183. https://doi.org/10.1007/s12028-011-9643-3